

การเลือกใช้งาน LC คอลัมน์และการแก้ปัญหาเบื้องต้น

จันทร์เพ็ญ ชัยวงศ์ขจร

บ.สิทธิพร แอสโซซิเอต จำกัด

21 กุมภาพันธ์ 2566

SITHIPORN
associates

HOW HPLC COLUMNS WORK

Our Products in Action!

phenomenex
...breaking with tradition

**LET'S ANALYZE
THE COMPOUNDS
IN COUGH AND
COLD REMEDIES.**



COUGH & COLD REMEDIES

2

SITHIPORN
associates


LC COLUMN TYPE

- NORMAL PHASE COLUMN
- REVERSE PHASE COLUMN
- ION EXCHANG COLUMN
- SIZE EXCLUSION COLUMN
- CHIRAL COLUMN


มีให้เลือกใช้มากที่สุด

REVERSE PHASE CHROMATOGRAPHY

Stationary phase = Non-Polar (C18)



- RED คือ สารที่ Polar มากที่สุด
- BLUE คือ สารที่ Polar รองลงมา
- GREEN คือ สารที่ polar น้อยที่สุด



COLUMN CHARACTERISTICS


- STRUCTURE OF SILICA
- BONDING TYPE
- PORE SIZE
- SURFACE AREA
- END CAPPING
- PARTICLE SIZE
- ID & COLUMN LENGTH

SILICA BASE

Fully Porous



Core-Shell



Two-In One (TWIN) Technology

การพัฒนาคุณภาพ ของ Siliga Base เพื่อ


- Pressure resistant
- pH Resistant
- Strong condition use
- Retention

SITHIPORN
associates

PORE SIZE

The smaller the pore size, the greater the surface area

การแยกของ Reverse Phase จะเกิดที่ Surface area



SITHIPORN
associates


BONDING CHEMISTRY

- C4, C6
- C8 >> USP L7
- C18 >> USP-L1-L2
- PHENYL >> USP L1 L
- CN >> USP L10
- NH2 >> USP L8
- อื่นๆ

USP Column Classification	Recommended Phases	Column	Particle Shape	Page
L1 Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod.	Kinetex C18		Core-Shell	235
	Kinetex EVO C18		Core-Shell	235
	Kinetex Polar C18		Core-Shell	235
	Kinetex PS C18		Core-Shell	235
	Kinetex XS-C18		Core-Shell	235
	Luna [®] C18Q2		Spherical	265
	Luna Omega C18		Spherical	279
	Luna Omega PS C18		Spherical	279
	Luna Omega Polar C18		Spherical	279
	Gemini [®] NX-C18		Spherical	223
	Gemini C18		Spherical	223
	Synergi [®] Hydro-RP		Spherical	332
	Synergi Fusion-RP		Spherical	332
	bioZen [®] Peptide-PS-C18		Spherical	208
	Onyx [®] C18		Monolith	302
	Juprelle [®] C18		Spherical	233
	Clarify [®] Oligo-RP		Spherical	394
	Clarify Oligo-MS		Core-Shell	394
	Clarify Oligo-TT		Core-Shell	394
	Aeris [™] WIDEPORE XB-C18		Core-Shell	202
	bioZen Peptide XB-C18		Core-Shell	208
L2 Octadecyl silane chemically bonded to silica gel of a controlled surface porosity that has been bonded to a solid spherical core, 30 to 50 µm in diameter.	Kinetex HILIC		Core-Shell	235
L3 Porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	Luna SilicaQ2		Spherical	265
	Onyx Silica		Monolith	302
L4 Silica gel of controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter.				
L5 Alumina of controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter.				
L6 Strong cation-exchange packing sulfonated fluorocarbon polymer coated on a solid spherical core, 30 to 50 µm in diameter.				
L7 Octyl silane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	Kinetex C8		Core-Shell	235
	Luna C8(2)		Spherical	265
L8 An essentially monomolecular layer of aminopropyl-silane chemically bonded to totally porous silica gel support, 1.5 to 10 µm in diameter, or a monolithic silica rod.	Onyx C8		Monolith	302
	bioZen Inertsil XB-C8		Spherical	208
L9 Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter.	Luna NH ₂		Spherical	265
	Luna Omega SUGAR		Spherical	279
	Luna SOX		Spherical	265
L10 Nitrile groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	Luna CN		Spherical	265
L11 Phenyl groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod.	Kinetex Esquire [®]		Core-Shell	235
	Kinetex Phenyl-Hexyl		Core-Shell	235

SITHIPORN
associates

ทบทวนความรู้พื้นฐาน VAN DEEMTER THEORY




SITHIPORN
associates

การเลือกคอลัมน์ มีปัจจัยอย่างไรบ้าง เมื่อ **BOND PHASE** เหมือนกัน

- PARTICAL SIZE
- ID COLUMN
- COLUMN LENGTH

SITHIPORN
associates

The Impact of Selectivity on Resolution




SITHIPORN
associates

COLUMN SIZE VS. EFFICIENTCY

Column Length (L) = 250 mm = 250,000 μ m

D_p = 5 μ m

$$N = \frac{L}{D_p} = \frac{250,000}{5} = 50,000$$



Length (L, mm)	Column Diameter (d_c , mm)	Particle Size (d_p , μ m)	L/d_p
250	4.6	10	25,000
150	4.6	5	30,000
150	2.1	5	30,000
100	4.6	3.5	28,600
100	2.1	3.5	28,600
75	4.6	2.5	30,000
75	2.1	2.5	30,000
50	4.6	1.7	29,400
50	2.1	1.7	29,400


SITHIPORN
associates

UPDATE COLUMN ใหม่ ๆ มีอะไรบ้างที่ต้องพิจารณา

Core-shell particles consist of a **solid core** coated with a layer of porous silica that is deposited either in layers or a single coating, depending on the manufacturer.


ประโยชน์ที่ได้ เมื่อเทียบกับ **Fully Porous** ทั่วไป
narrower particle size distribution
significantly lower back pressures

ปัจจุบัน **Core shell column** มีหลากหลายให้เลือกทุกยี่ห้อ มีตั้งแต่ **Analytical** ถึง **Preperative scale**




Retention on Core-Shell Media

Surface Area




UHPLC Performance – Cyclosporine Impurity Profile


Kinetex 2.6µm Polar C18



Luna Omega 1.6µm Polar C18



Conventional 1.7µm C18



Conditions for all columns same except where noted:

Columns: Kinetex 2.6µm Polar C18

Luna Omega 1.6µm Polar C18

Conventional Fully Porous 1.7µm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: Acetonitrile/Tert-butyl methyl ether/Water/
Phosphoric acid (430:50:520:1)

Flow Rate: 0.30 mL/min

Temperature: 80 °C


Detection: UV @ 210 nm

Sample: Cyclosporine


SITHIPORN
associates

DECREASE PARTICLE SIZE- LOWER LEVEL OF DETECTION AND QUANTITATION

Agilent Technologies® ZORBAX® 3.5µm SB-C18



Kinetex 2.6µm C18



Conditions for both columns:

Dimensions: 150 x 4.6 mm

Mobile Phase: A: Water

B: Acetonitrile

Gradient: (95:5) A/B for 1.16 min, then to (5:95) A/B

Flow Rate: 1.5 mL/min

Temperature: 45 °C







Detection: UV @ 254 nm

Instrument: Agilent 1200

Sample:


- | | |
|---------------------|-------------------------------------|
| 1. Pyridine | 9. Nortriptyline |
| 2. Acetaminophen | 10. 4-Chlorobenzoic acid |
| 3. Pindolol | 11. 5-Methyl-2-hydroxy benzaldehyde |
| 4. Quinine | 12. 4-Chlorocinnamic acid |
| 5. Acebutolol | 13. Diazepam |
| 6. Chlorpheniramine | 14. Diflunisal |
| 7. Triprolidine | 15. Niflumic acid |
| 8. Prednisolone | 16. Hexanophenone |

SITHIPORN
associates

					
100% Aqueous , good peak shape for basic compound	Polar & Non-Polar retention , for 100% Aqueous mobile phase	High pH stable C18	Improved peak shape of basic compound and increased retention of acidic compound	Core-Shell C18 universal	Core-Shell C8 universal
pH Range 1.5-8.5 USP L1	pH Range 1.5-8.5 USP L1	pH Range 1-12 USP L1	pH Range 1.5-8.5 USP L1	pH Range 1.5-8.5 USP L1	pH Range 1.5-8.5 USP L7

RELATING SELECTIVITY TO UHPLC STATIONARY PHASE

Define the Characteristics of Your Target Compounds



Hydrogen Bond Accepting Capacity

Hydrophobicity


π - π Interaction

SITHIPORN
associates

POLAR SELECTIVITY OF CATECHOLAMINES


Conditions for both columns:

Columns: Luna Omega 1.6µm Polar C18
Kinetex 1.7µm C18
Dimensions: 50 x 2.1 mm
Mobile Phase: A: 10 mM Ammonium Formate with 0.1% Formic Acid
B: Acetonitrile with 0.1% Formic Acid
Gradient: Time (min) % B
0 0
3 90
Flow Rate: 0.4 mL/min
Injection Volume: 1 µL
Temperature: 22°C
Detection: MS/MS (SCIEX API 4000™)
Sample: 1. Norepinephrine
2. Epinephrine
3. Normetanephrine
4. Dopamine
5. Metanephrine
6. Serotonin




SITHIPORN
associates


Kinetex 5µm EVO C18



Waters® XBridge® 3.5µm C18



Waters XBridge 5µm C18



Conditions for both columns:

Column: Kinetex 5µm EVO C18
XBridge 5µm C18
XBridge 3.5µm C18
Dimensions: 150 x 4.6 mm
Mobile Phase: A: 0.1% Formic acid in Water
B: 0.1% Formic acid in Acetonitrile
Gradient: 5% to 95% B over 10 minutes
Flow Rate: 1.5 mL/min
Temperature: Ambient
Detection: UV @ 254 nm

Sample: 1. Uracil
2. Pindolol
3. Chlorpheniramine
4. Nortriptyline
5. 3-Methyl-4-Nitrobenzoic acid
6. 5-Methyl Salicyl Aldehyde
7. Hexaphenone

SITHIPORN
associates

Avoiding Problems

“Chromatographic Hygiene”

- GUARD COLUMNS
 - Particulate matter from piston seals and injection valve rotors can clog columns if they are not removed
- 0.2 μm SYRINGE FILTER
 - With sub-2 μm columns, every sample must be filtered through a 0.2 μm porosity filter
- MOBILE PHASE FILTERS
 - Fresh buffers must be made daily and filtration through a 0.2 μm porosity filter is required
- INLINE FILTERS:
 - I strongly recommend using a 0.2 μm porosity in-line filter between the autosampler and guard column

Benefits of Guard Columns

- Why Use Guard Columns?
 - Protect valuable analytical columns by removing particulates and strongly retained sample components that may accumulate on column
 - Increase lifetime of analytical column
 - Maintain
 - High column efficiencies
 - Resolution
 - Peak shape
 - Cost-effective


Benefits of Guard Columns

Protected Column

Inlet Frit



Column Media




Unprotected Column

Enrichment of particles and impurities

Inlet Frit



Column Media




 *...breaking with tradition™*

SecurityGuard™

SecurityGuard™ ULTRA





- SecurityGuard (left) for fully porous $\geq 3 \mu\text{m}$ columns
 - Cartridge only 4 mm long
- SecurityGuard ULTRA (right) for core-shell and sub-2 μm columns (2.0 to 4.6 mm ID)
 - Low dead volume ($<0.3 \mu\text{L}$)
 - Pressure rated up to 20,000 psi (1378 bar)
- Protects against damaging chemical contaminants and microparticulates
- Will NOT alter chromatography & Easy to Use

 *...breaking with tradition™*

When Should I Replace Guard Column?

Symptoms	Criteria
• Increasing system backpressure	• $>20 \%$ or before system auto-shutdown
• Loss of peak efficiency	• Efficiency (N) decreases $>20 \%$
• Loss of peak resolution (merging, shifting)	• Resolution (R_s) decreases 10%
Every time the analytical column is replaced, the guard column should also be replaced.	





phenomenex
...breaking with tradition™


Mobile Phase Filtration

Phenex™ Filter Membrane

- 47 mm diameter
- 0.2 + 0.45 μm
- RC, PTFE, PES, NY

Funnel for successive filtration



Solvent pickup adaptor for continuous filtration

SITHIPORN associates

phenomenex
...breaking with tradition™

Solvent Inlet Filter

- Unfiltered or contaminated solvents (microbiological growth) can clog solvent inlet filters
- This has an impacting on pump performance
- HPLC system has to work harder to pull mobile phase into the pump or even sucks air
- Clean the solvent inlet filters whenever preparing new mobile phase and replace it once every 3-6 months



SITHIPORN associates


phenomenex
...breaking with tradition™

Mobile Phase Filtration

- Strongly recommend preparing aqueous buffers fresh DAILY to minimize microbial growth, which can foul in-line filters and columns

Microbial growth on packing material at inlet of column

SITHIPORN
associates



phenomenex
...breaking with tradition™


Sample Filtration

Syringe Filter




Luer lock


Membrane sealed within a polymeric housing

Outlet



Select filter diameter based on sample volume

If your sample volume is:		
< 2 mL Sample Volume	2 to 10 mL Sample Volume	10 to 100 mL Sample Volume
4 mm Diameter	15 - 17 mm Diameter	25 - 28 mm Diameter
		



Particulates Trapped

Particulate-free Sample

Preparing Isocratic Mobile Phase

Isocratic Mobile Phase Preparation

- What is the correct way to mix isocratic mobile phase; and how does it impact chromatography?
- Let's run the following experiment and see if it matters how an isocratic mobile phase is prepared:
 - A. Pre-mix
 - i. Measure water volume, then add organic to total volume?
 - ii. Measure organic volume, then add water to total volume?
 - iii. Measure volume for each component separately, then mix?
 - B. Allow HPLC pump to mix A and B and deliver to column



Isocratic Mobile Phase Preparation


- Pre-Mix
 - Water first:
 - 400 mL water in graduated cylinder, add to 1-Liter volumetric flask, dilute to mark with MeOH
 - MeOH first:
 - 600 mL MeOH in graduated cylinder, add to 1-Liter volumetric flask, dilute to mark with Water
 - MeOH and Water measured individually:
 - 600 mL MeOH in graduated cylinder, transfer mobile phase reservoir; 400 mL Water in same graduated cylinder and transferred to mobile phase reservoir
- Allow pump to mix A = Water with B = MeOH

SITHIPORN
associates



Isocratic Mobile Phase Preparation

Kinetex 5 μ m XB-C18
150 x 4.6
60:40 MeOH : H₂O
Flow rate: 1.25 mL/min
Inj. Vol: 2 μ L
UV @ 254 nm



SITHIPORN
associates

phenomenex[®]
...breaking with tradition[™]

Troubleshooting Approach & Diagnostic Tools

SITHIPORN
associates


phenomenex[®]
...breaking with tradition[™]

Tailing – All Peaks


Diagnostic Considerations

All Peaks = Mechanical

1. Column Void (gap)
2. Connectivity
3. Physical/Particulate Blockage




SITHIPORN
associates



phenomenex
...breaking with tradition™

Column Frits

1. Contain media within column while allowing mobile phase to pass
2. Larger particles can block flow path
3. Peak Tailing and High Backpressure





SITHIPORN
associates

Tailing – Isolated Peak

Diagnostic Considerations

Isolated Peak(s) = Chemical Interaction

1. Free (weak acid) silanols
2. Cation-Exchange with Basic Compounds
3. Secondary Interaction Mechanism



Column Care



Column Cleaning

- I am often asked “how can I clean my HPLC/UHPLC column?”
- My typical answer is: “do not get it dirty”
- A bit flippant, however, so let me explain.....

SITHIPORN
associates



Column Cleaning

- Consult with column manufacturer for recommended cleaning procedures and limitations
 - Reverse flushing is best if allowed
- Flush column with mobile phase minus buffer to remove all buffer salts
- Flush with 10 column volumes (or more) of mutually miscible solvent such as methanol or acetonitrile
- Flush with 20 column volumes (or more) of strong solvent (THF or IPA) to remove strongly adsorbed sample components
- Reverse the process and equilibrate with mobile phase & retest the column

SITHIPORN
associates



Routine Column Cleaning

Reversed phase silica columns

1. Flush with water:organic (methanol or acetonitrile) 95:5 or the same ratio as used in the analysis
2. Increase the organic percentage in a gradient up to 100% (or at least 20% greater than final mobile phase ratio)
3. Flush with 40:60 water:organic or any similar ratio for short-term storage

SITHIPORN
associates



Column Storage

Reversed phase silica columns

1. Remove buffer salts before storage
2. Keep pores wetted with at least 30% organic
3. Elute adsorbed components before storage, but do not precipitate out buffer salts in the column
4. Never leave a column static at an elevated temperature

SITHIPORN
associates



Tips for Extending Column Lifetimes & Minimizing Problems

- **Make sure system is clean before installing a new column**
 - This includes flushing old mobile phase and buffer from system before installing new column
- **Make fresh mobile phase daily**
 - Microbacterial growth can occur in aqueous mobile phase, and this can damage columns and systems
- **Use guard columns**
- **Filter samples before injection, no matter how 'clean' they appear**

SITHIPORN
associates



Summary

- A. Practice Good Chromatographic Hygiene**
- B. Check the Simple Things First**
- C. Structured Diagnostic Approach**
 - 1. Mechanical pathway
 - 2. Chemical interaction pathway
 - 3. Diagnostic tools – Full Chromatographic Profile, Backpressure, Void Time (t_0)
- D. Column Care**

SITHIPORN
associates