




Nitrosamines: Overview & Solutions

Ajiraporn Kongpol, Ph.D Country Sales Manager

22 March 2024



Nitrosamine Impurities in Drugs



- Nitrosamines, or N-nitrosamines, are chemical compounds believed to be carcinogenic to humans based on animal studies.
- Therefore, nitrosamines need to be highly regulated to ensure that there is no adverse effect to human health.
- In 2018, multiple nitrosamines including NDMA were detected within various blood pressure medicines known as ‘sartans’.
- Since then, nitrosamines have been detected in ranitidine, metformin, and several other drugs which highlighted that the issue was not solely limited to sartans.
- The issue has resulted in multiple product recalls worldwide, with the problem likely to continue for some time as more drug products become implicated.
- Regulatory agencies (FDA, EMA, USP, Health Canada, etc.) were quick to act in establishing guidance for nitrosamine limits, testing, and risk assessment for prevention of nitrosamine formation in drug substance (API) and drug products.

Nitrosamine Formation in Drugs

Figure 1. Representative Reaction to Form Nitrosamines





Figure 3. Formation of NDMA From N,N-Dimethylformamide



- Formation of nitrosamines is possible in the presence of secondary, tertiary, or quaternary amines and nitrite salts under acidic reaction conditions. Under these conditions, nitrite salts may form nitrous acid, which can react with an amine to form a nitrosamine (see Figure 1).
- Amide solvents, which are susceptible to degradation under certain reaction conditions, are another source of secondary amines. For example, under high reaction temperatures for an extended reaction period, N,N-dimethylformamide can degrade into dimethylamine, which can react with nitrous acid to form NDMA (see Figure 3).

USP General Chapter <1469>

- 7 Nitrosamine Reference Standards available
- 4 Analytical Procedures
 - Quantitation of NDMA, NDEA, NDIPA, NEIPA, NMBA, NMPA, and NDBA in selected sartans (valsartan, irbesartan, and losartan potassium) by HPLC–HRMS on **Kinetex F5 (L43)** 2.6 µm 100 x 4.6 mm.
 - Quantitation of NDMA, NDEA, NDIPA, and NEIPA in selected sartans (valsartan, irbesartan, losartan potassium, olmesartan medoxomil, candesartan cilexetil, and telmisartan) by headspace GC–MS on **Stabilwax (G16)** 30 m x 0.32 mm x 1.0 µm.
 - Quantitation of NDMA, NDEA, NDIPA, NEIPA, NMBA, and NDBA in selected sartans (valsartan, losartan potassium, olmesartan medoxomil, candesartan cilexetil, and telmisartan) by HPLC–MS/MS on **Raptor ARC-18 (L1)** 2.7 µm 150 x 3.0 mm.
 - Quantitation of NDMA, NDEA, NDIPA, NEIPA, NMPA, and NDBA in selected sartans (valsartan, losartan potassium, and candesartan cilexetil) by GC–MS/MS (triple-quad) on **VF-WAXms (G16)** 30 m x 0.25 mm x 1.0 µm.

NOTE: ALL published methods use MS for the low-level detection required

FDA Methods

- FDA has published several methods, including 2 on PHX columns



Liquid Chromatography-High Resolution Mass Spectrometry (LC-HRMS) Method for the Determination of Six Nitrosamine Impurities in ARB Drugs

HPLC column: Kinetex® 2.6 µm F5 100 Å, 100 x 4.6 mm (Phenomenex, Part No. 00D-4723-E0)

06/03/2020



HPLC column: Phenomenex Kinetex® 2.6 µm Biphenyl 100 Å, 150 x 3.0 mm (Part No. 00F-4622-Y0)

Liquid Chromatography-Electrospray Ionization-High Resolution Mass Spectrometry (LC-ESI-HRMS) Method for the Determination of Nitrosamine Impurities in Metformin Drug Substance and Drug Product

NOTE: ALL published methods use MS for the low-level detection required

Nitrosamine Analytical Challenges

- Nitrosamines lack a chromophore
- Low-level (ppb) detection and quant of nitrosamines is required
 - UV detectors (workhorse of Pharma QC labs) are NOT sensitive enough
 - MS detection required
- High concentration of API in drug substance and drug product
- No single or generic method is available
- Chromatographic interference
 - Separation of API from individual nitrosamines is required in drug product and drug substance to enable accurate quantitation


Analyzing Nitrosamines

Nitrosamines are known to be carcinogenic and have the potential to be intermediates in organic synthesis. Due to their potent genotoxicity, nitrosamines have been a serious cause for concern. It has now become a requirement to accurately quantitate this group of compounds in pharmaceuticals during drug development and manufacturing. Below you will find the latest resources to help you improve your LC-MS quantitation of nitrosamines. If you are looking to improve your nitrosamine analysis please do not hesitate to contact us.


| Technique | Brand | Description | Recommended for: |
|-----------|------------|--|---|
| HPLC | Gemini | 3 μ m C18 150 x 4.6 mm | Ph. Eur. method for N-Nitrosamines in active substances |
| HPLC | Kinetex | 2.6 μ m Biphenyl 150 x 3 mm 2.6 μ m Biphenyl 150 x 4.6 mm | Sartan drugs, metformin and ranitidine |
| HPLC | Kinetex | 2.6 μ m F5 50 x 2.1 mm | Rapid analysis of genotoxic nitrosamines |
| GC | Zebron | ZB-624PLUS 30 m x 0.25 mm x 1.40 μ m | Analysis of NDMA and NDEA in sartan drugs |
| HPLC | Luna Omega | 3 μ m Polar C18 150 x 3 mm | Nitrosamine impurities in acyclovir |

PLC/UHPLC Column Phase Selection Chart

Specifically designed for successful and reproducible method development and transfer



Gemini Twin Technology



Gemini-NXTwin Technology

Gemini columns are rugged reversed-phase HPLC and PREP LC columns that offer extended lifetime at pH 1-12 conditions and excellent stability for reproducible, high-efficiency separations. The Gemini family contains a variety of useful alkyl and phenyl phases bonded to fully porous organo-silica particles.



Rugged reversed phase column can be used reliably from pH 1-12

| Phase | | Particle Type | Pore Size (Å) | Carbon Load % | pH Range | USP Classification | Selectivity Features | Available Particle Size(s) |
|------------------|--|--|---------------|---------------|----------|--------------------|--|----------------------------|
| C18 | <p>Hydrophobicity</p> <p>Steric Interaction</p> <p>Hydrogen Bond Donating Capacity</p> <p>Hydrogen Bond Accepting Capacity</p> <p>Cation Selectivity at pH 2.8</p> <p>Cation Selectivity at pH 7.0</p> <p>Low High</p> | Fully Porous Organo-Silica (U.S. Patent Nos. 7,563,367 and 8,658,038) | 110 | 14 | 1.0 - 12 | L1 | This is a high loading, organo-silane particle column with pH stability 1-12. The patented procedure creates a surface that is a strong hydrogen donor and acceptor. It is ideal for a combination of polar and non-polar retention. | 3µm 5µm 10µm |
| NX-C18 | <p>Hydrophobicity</p> <p>Steric Interaction</p> <p>Hydrogen Bond Donating Capacity</p> <p>Hydrogen Bond Accepting Capacity</p> <p>Cation Selectivity at pH 2.8</p> <p>Cation Selectivity at pH 7.0</p> <p>Low High</p> | Fully Porous Organo-Silica (U.S. Patent Nos. 7,563,367 and 8,658,038) | 110 | 14 | 1.0 - 12 | L1 | New generation of organo-silane material incorporates ethylene bridges to provide pH stability from 1-12 and 5x the durability of earlier hybrids. The homogenous surface offers some steric selectivity. | 3µm 5µm 10µm |
| C6-Phenyl | <p>Hydrophobicity</p> <p>Steric Interaction</p> <p>Hydrogen Bond Donating Capacity</p> <p>Hydrogen Bond Accepting Capacity</p> <p>Cation Selectivity at pH 2.8</p> <p>Cation Selectivity at pH 7.0</p> <p>Low High</p> | Fully Porous Organo-Silica (U.S. Patent Nos. 7,563,367 and 8,658,038) | 110 | 12 | 1.0 - 12 | L11 | This is a very inert phase for great peak shapes of ionized compounds. The planar phenyl rings offer moderate hydrophobic retention and high steric selectivity for structural isomer selectivity. | 3µm 5µm |

PLC/UHPLC Column Phase Selection Chart




Use for highest performance compared to fully porous media

| Phase | | Particle Type | Pore Size (Å) | Carbon Load % | pH Range | USP Classification | Selectivity Features | Available Particle Size(s) |
|--|--|--|---------------|-----------------|----------------------|--------------------|---|----------------------------|
| EVO C18 | | Core-Shell Organo-Silica (U.S. Patent Nos. 7,563,367 and 8,658,038) | 100 | 11 ^a | 1.0-12.0 | L1 | Robust reversed phase methods even in alkaline conditions with improved peak shape for polar basic compounds and 100% aqueous stability. | 1.7 µm 2.6 µm 5 µm |
| C18 | | Core-Shell Silica | 100 | 12 ^a | 1.5-8.5 ^a | L1 | All purpose phase that offers the hydrophobic retention and methylene selectivity chromatographers expect from a C18 column. | 1.3 µm 1.7 µm 2.6 µm 5 µm |
| Biphenyl | | Core-Shell Silica | 100 | 11 ^a | 1.5-8.5 ^a | L11 | 100% aqueous stable and allows for excellent reversed phase retention and enhanced polar and aromatic selectivity. | 1.7 µm 2.6 µm 5 µm |
| XB-C18 (C18 with protective iso-butyl side chains) | | Core-Shell Silica | 100 | 10 ^a | 1.5-8.5 ^a | L1 | Unique C18 phase that yields increased hydrogen bonding with hydrophobic selectivity, resulting in improved peak shape for basic compounds and increased retention of acidic compounds. | 1.7 µm 2.6 µm 3.5 µm 5 µm |
| F5 (pentafluorophenyl propyl) | | Core-Shell Silica | 100 | 9 ^a | 1.5-8.5 ^a | L43 | Highly reproducible pentafluorophenyl propyl phase that offers a unique combination of polar, hydrophobic, aromatic, and shape selectivity. | 1.7 µm 2.6 µm 5 µm |
| C8 | | Core-Shell Silica | 100 | 8 ^a | 1.5-8.5 ^a | L7 | USP L7 phase that provides less hydrophobic and methylene selectivity than a C18. | 1.7 µm 2.6 µm 5 µm |
| Phenyl-Hexyl | | Core-Shell Silica | 100 | 11 ^a | 1.5-8.5 ^a | L11 | Reversed phase chemistry that allows for greater retention and separation of aromatic hydrocarbons. | 1.7 µm 2.6 µm 5 µm |
| HILIC (Unbonded silica) | | Core-Shell Silica | 100 | 0 | 2.0-7.5 | L3 | Unbonded silica phase for HILIC conditions to provide selectivity for polar compounds. | 1.7 µm 2.6 µm 5 µm |

^a Effective Carbon Load. ^a pH range is 1.5-10 under isocratic conditions. pH range is 1.5-8.5 under gradient conditions.



Kinetex Core-Shell Technology delivers dramatic improvements in efficiency over conventional fully porous media, which can be leveraged to increase resolution, greatly improve productivity, reduce solvent consumption, and decrease costs.



Using sol-gel processing techniques that incorporate nano-structuring technology, a durable, homogeneous porous shell is grown on a solid silica core to create a Kinetex Core-Shell particle. This particle morphology results in less band broadening for all four sources of UHPLC band broadening compared to fully porous particles and thus delivers extremely high efficiencies.




One of the world's leading HPLC/UHPLC columns for virtually every application




Luna Omega Silica

One of the world's leading HPLC brands, now enhanced for incredible HPLC and UHPLC performance! Luna Omega 1.6, 3, and 5 μm columns culminate 20 years of technological prowess, advancements, and innovation from Phenomenex! With astounding efficiency levels, highly versatile selectivities, and trusted accuracy, Luna Omega columns will take your chromatographic experience to a new level.



Through our proprietary process, we eliminate micropores, further improving column efficiency, inertness, and reproducibility.



Consistent Porosity

Absence of Micropores

| Phase | Particle Type | Pore Size (\AA) | Carbon Load % | pH Range | USP Classification | Selectivity Features | Available Particle Size(s) |
|---|--|----------------------------|-------------------|----------------------|--------------------|--|--|
| Omega C18 Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Thermally Modified Fully Porous Silica | 100 | 11 | 1.5-8.5 ^a | L1 | Rugged and highly efficient C18 with strong focus on hydrophobic retention of non-polar and polar compounds | 1.6 μm |
| Omega Polar C18 Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Thermally Modified Fully Porous Silica | 100 | 9 | 1.5-8.5 ^a | L1 | 100% aqueous stability and enhanced selectivity/retention for polar analytes without diminishing useful non-polar retention. The C18 ligand provides general hydrophobic interactions while a polar modified particle surface provides enhanced polar compound retention. | 1.6 μm 5 μm |
| Omega PS C18 Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Thermally Modified Fully Porous Silica | 100 | 9 | 1.5-8.5 ^a | L1 | Unique, 100% aqueous stable mixed-mode phase that provides both polar and non-polar retention. The surface contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while the C18 ligand promotes general reversed phase retention. The positively charged surface also improves basic compound peaks shape through ionic repulsion. | 1.6 μm 5 μm |
| C18(2) Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 17.5 | 1.5-9.0 ^a | L1 | Excellent general purpose reversed phase selectivity with high hydrophobic and methylretent. Non-polar endcapping virtually eliminates silanol interactions. | 2.5 μm 3 μm 5 μm 10 μm |
| C8(2) Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 13.5 | 1.5-9.0 ^a | L7 | General C8 phase that provides less hydrophobic retention than a C18, but the density of the ligand bonding creates more steric based selectivity. | 3 μm 5 μm 10 μm |
| CN Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 7 | 1.5-7.0 | L10 | Can be used as reversed or normal phase material. Nitrile groups bound to the silica surface offer a unique polar selectivity under reversed phase or normal phase conditions. | 3 μm 5 μm 10 μm |
| HILIC (cross-linked diol) Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 200 | 5.7 | 1.5-8.0 | L20 | HILIC phase that provides excellent selectivity and retention for polar compounds as well as improved MS sensitivity with low bleed. | 3 μm 5 μm |
| NH₂ Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 9.5 | 1.5 - 11 | L8 | Amino phase. Can be used in reversed or normal phase modes. Stable from pH 1.5 to 11.0 and under 100% aqueous conditions. High performance silica and bonding techniques produces a rugged, highly reproducible column. | 3 μm 5 μm 10 μm |
| PFP(2) Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 11.5 | 1.5-8.0 | L43 | A pentafluorophenyl phase that provides excellent selectivity for aromatic compounds from influence of fluorine substitution on phenyl ring. Multiple retention mechanisms. Orthogonal selectivity to traditional C18 phases. | 3 μm 5 μm |
| Phenyl-Hexyl Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 17.5 | 1.5-9.0 ^a | L11 | A phenyl phase which employs a hexyl alkyl linker, leading to a balanced degree of hydrophobic and aromatic selectivity. Our most hydrophobic phenyl column and it will also provide good hydrogen accepting functionality for acidic retention. | 3 μm 5 μm 10 μm |
| SCX Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 0.55% Sulfur Load | 2.0 - 7.0 | L9 | A benzoate sulfonic acid bonded phase is used to make this strong cation-exchange (SCX) column. Offers great peak shape and resolution. | 5 μm 10 μm |
| Silica(2) Hydrophobicity Steric Interaction Hydrogen Bond Donating Capacity Hydrogen Bond Accepting Capacity Cation Selectivity at pH 2.0 Cation Selectivity at pH 7.0 | Fully Porous Silica | 100 | 0 | 2.0 - 7.5 | L3 | Ultra-pure silica with high column bed stability enhanced by particle shape uniformity. | 3 μm 5 μm 10 μm |

TN-1334

Sensitive and Reproducible Quantification of N-Nitroso Propranolol in a Propranolol Drug Substance and Product

Pankaj Partani¹, Sandeep Choudhary, PhD¹, Preeti Bharatiya², Upendra Gunta³, Ranjith Kumar Ponnannaneni³, Manoj Pillai, PhD¹, Rahul Baghla², Eshani Nandita, PhD², and Bryan Tackett, PhD³


¹SCIEX Lab, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Hobli, Jala Taluka, Bengaluru 562149 India

²AB Sciex LLC, 500 Old Connecticut Path, Framingham, MA 01701, USA


³Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA



1.0 CC Std#07#19 - N-Nitroso Propranolol...ity\Sample.078.wiff), (sample Index: 1)
Area: 864915, Height: 160413, RT: 7.322 min



Rec-Tab-1.0#01 - N-Nitroso Propranolol...ity\Sample.079.wiff), (sample Index: 1)
Area: 1010077, Height: 188281, RT: 7.319 min



LC Conditions

Column: Kinetex™ 2.6 µm Biphenyl

Dimensions: 150 x 3.0 mm

Part No.: [00F-4622-Y0](#)

Mobile Phase: A: 1 mM Ammonium Formate with 0.1 % Formic Acid in Water
B: 0.1 % Formic Acid in Acetonitrile

| Gradient | Time (min) | %B |
|----------|------------|----|
| | 0.00 | 35 |
| | 3.00 | 35 |
| | 5.00 | 55 |
| | 6.00 | 75 |
| | 8.00 | 75 |
| | 8.10 | 95 |
| | 10.0 | 95 |
| | 10.1 | 35 |
| | 14.0 | 35 |

Flow Rate: 0.4 mL/min

Injection Volume: 15 µL

Temperature: 40 °C

LC System: SCIEX ExionLC™

Detection: MRM

Detector: SCIEX QTRAP 6500+

TN-1335



Low-level Quantification of 10 Mutagenic Nitrosamine Impurities in Acyclovir

Lakshmanan Deenadayalan¹, Sashank Pillai¹, Rahul Baghla², Elliot Jones², Eshani Nandita, PhD², and Bryan Tackett, PhD³

¹SCIEX Lab, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Hobli, Jala Taluka, Bengaluru 562149

²AB Sciex LLC, 500 Old Connecticut Path, Framingham, MA 01701, USA

³Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA



LC Conditions

Column: Luna™ Omega 3.0 μm Polar C18

Dimensions: 150 x 3.0 mm

Part No.: [00F-4760-Y0](#)

Mobile Phase: A: 0.1 % Formic Acid in Water

B: 0.1 % Formic Acid in Methanol

| Gradient: | Time (min) | %B |
|-----------|------------|----|
| | 0.00 | 2 |
| | 1.00 | 2 |
| | 20.0 | 98 |
| | 22.0 | 98 |
| | 22.1 | 2 |
| | 24.0 | 2 |

Flow Rate: 0.45 mL/min

Injection Volume: 20 μL

Temperature: 35 °C

LC System: SCIEX® ExionLC™

Detection: MRM

Detector: SCIEX QTRAP® 6500+

LC-MRM-MS Method for the Detection and Quantification of Six Nitrosamine Impurities in Sartan (ARBs) Drugs

Results

Figure 1. Representative Chromatogram of 6 Nitrosamine Impurities Using a Kinetex 2.6 μ m F5 Column.




Figure 4. Calibration Curve for the Six Nitrosamines Impurities Using a Kinetex 2.6 μ m F5 Column. Calibration Range is from 0.5 to 25 ng/mL with Correlation Coefficient >0.980 and LOD of 0.007 ppm.




Figure 2. Representative Chromatogram of 6 Nitrosamine Impurities Using a Kinetex 2.6 μ m Biphenyl Column.




Figure 5. Calibration Curve for the Six Nitrosamines Impurities Using a Kinetex 2.6 μ m Biphenyl Column. Calibration Range is from 0.5 to 25 ng/mL with Correlation Coefficient >0.998 and LOD of 0.005 ppm.




Figure 3. Representative Chromatogram of 6 Nitrosamine Impurities Using a Luna Omega 3 μ m Polar C18 Column.





Figure 6. Calibration Curve for the Six Nitrosamines Impurities Using a Luna Omega 3 μ m Polar C18 Column. Calibration Range is from 0.5 to 25 ng/mL with Correlation Coefficient >0.980 and LOD of 0.005 ppm.



LC-MRM-MS Method for the Detection and Quantification of Six Nitrosamine Impurities in Sartan (ARBs) Drugs

Results

Figure 7. % Recovery of Nitrosamine Impurities Spiked in 100 mg/mL Valsartan Using a Kinetex 2.6 μ m F5 Column. Recovery of three replicates fell between 100 and 202 % with %RSD between 1.40 and 12.5 %.




Figure 8. % Recovery of Nitrosamine Impurities Spiked in 100 mg/mL Valsartan Using a Kinetex 2.6 μ m Biphenyl Column. Recovery of three replicates fell between 110 and 135 % with %RSD between 18.5 and 5.6 %.





Figure 9. % Recovery of Nitrosamine Impurities Spiked in 100 mg/mL Valsartan Using a Luna Omega 3 μ m Polar C18 Column. Recovery of three replicates fell between 100 and 202 % with %RSD between 1.40 and 12.5 %.



Discussion

Conventional C18 columns provide little interaction and retention for small, highly polar compounds. Nitrosamines are compounds of low molecular weight and high polarity, thus making chromatographic separation a challenge. It is important for the column to be capable of retaining the polar Nitrosamines without generating secondary interactions. Secondary interactions can affect peak shape with issues such as tailing or reduced peak height. In turn, this will impact the sensitivity and quantitation of the analytical method.

Conclusion

Here, we present LC-MRM-MS methods capable of detecting highly polar Nitrosamines at 0.007 ppm in the drug product, which is lower than the specification limit for NDMA. Three different stationary phases, a highly inert pentafluoro phenyl, an aqueous stable biphenyl, and a Polar C18 phase, capable of providing excellent retention of polar Nitrosamines with good peak shape, were used for this purpose.

Phenomenex & SCIEX Tech Notes

- [TN-1259: Rapid Analysis of Genotoxic Nitrosamines by HPLC-MS/MS](#)
- [TN-1265: LC-MS/MS Quantitative Analysis of NDMA in Ranitidine using Kinetex Biphenyl](#)
- [Quantification of genotoxic nitrosamines in a telmisartan drug product](#)
- [Rapid analysis of genotoxic nitrosamines by HPLC-MS/MS](#)
- [Nitrosamine analysis in a pioglitazone drug product](#)
- [Rapid method for quantifying nitrosamine compounds with qualitative confirmation](#)
- [Analysis of nitrosamine impurities in a metformin drug substance and drug product](#)
- [Precise and accurate quantification of nitrosamine impurities in an esomeprazole API](#)
- [Varenicline nitrosamine drug substance-related impurity \(NDSRI\) quantification in a varenicline drug product](#)
- [LC-MRM-MS Method for the Detection and Quantification of Six Nitrosamine Impurities in Sartan \(ARBs\) Drugs](#)
- [Robust Analysis of Nitrosamines in a Losartan Drug Substance using Strata™ Activated Carbon Extraction](#)

Extraordinary Technical Support

AVAILABLE 24/7

LIVE CHAT WITH TECHNICAL EXPERTS

- 24/7 globally

CONCIERGE-STYLE SUPPORT

- Personalized programs

EASY ACCESS TO INFORMATION

- Applications, webinars, tech notes, white papers
- New website

EASY ORDERING AND HIGH ON-TIME DELIVERY STANDARDS

- Standing orders, subscription ordering


SHOULDER-TO-SHOULDER [EVEN IN A VIRTUAL WORLD]

- We want to be part of your team



Local and Global Support

- Over 1,100 employees in 25 offices worldwide and growing
- Directly servicing over 100 countries
- 10 manufacturing lines in 4 plants located in California, Colorado, and China.
 - Including U/HPLC, GC, SPE (silica), SPE (polymer), Proficiency Standards, Bulk media, Flash, and Instrumentation
- 8,500+ products
- 11 R&D 100 awards for innovation
- 92% customer satisfaction



Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at <https://www.phenomenex.com/phx-terms-and-conditions-of-sale>.

Trademarks

Kinetex and Luna are trademarks of Phenomenex.

SCIEX and QTRAP are registered trademarks and ExionLC is a trademark of AB SCIEX Pte. Ltd.

Disclaimer

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.